

# NGS Bioinformatics

**Module topic:** NGS Data Processing and QC

**Contact session title:** Read Data QC

**Trainer:** Shaun Aron

**Participant:** <write your name here>

**Date:** <write today's date here>

## FASTQ Data QC

### Introduction

In the assignment, you will have an opportunity to view and run a set of FASTQ files through the FASTQC program and assess the html output reports.

- View raw high throughput sequencing files (FASTQ).
- Check the quality of the dataset = run a standard QC pipeline.
  - This step is crucial to ensure that your dataset is of good quality before you continue further with the analysis.
- Remove adapters and low quality reads from the dataset.
- Compare a « good » and a « bad » quality dataset.

### Tools used in this session

FASTQC - <https://www.bioinformatics.babraham.ac.uk/projects/fastqc/>

Trimmomatic - <https://github.com/usadellab/Trimmomatic>

### **Task 1: Quality Control of a FASTQ dataset**

*The datasets for this assignment can be found on the course website. The dataset you will be working with first is paired end reads from an Escherichia coli K1 sample. The data has been produced on the Illumina MiSeq platform and is stored in the SRR957824.zip archive.*

Log into the HPC

In your home directory create a folder called ngs\_qc

```
mkdir ngs_qc  
cd ngs_qc
```

Download the datasets into the ngs\_qc folder

```
wget https://hpc.ilri.cgiar.org/~mkofia/NGS_QC/SRR6319976.zip  
wget https://hpc.ilri.cgiar.org/~mkofia/NGS_QC/SRR957824.zip
```

Unzip the SRR957824 archive and cd into the unzipped directory

```
unzip SRR957824.zip
```

```
cd SRR957824
```

View the files in the directory.

```
ls
```

View the FASTQ files

```
zcat SRR957824_1.fastq.gz | less
```

Tip: You can use the *space bar* to scroll through the file and *Q* to quit. The suffix of *\_1* indicates that the file contains data for Read1.

```
@SRR957824.1 1/1
ACTATGAGCGAAACGGCATCCTGGCAGCCGAGCGCATCCATTCTAACTTATTAACGCGCGCGATTATGGCGGAGAACCTCGTTCTTTCCCATCGTGGAGTGCTGGAGGTGGAGACGCCCTGTATGTGCCAGGCGACGGTAAC
+
?????BBBDDDDDBDFFFFFFFHHHICEEEHBHHIIIIIFGHIIIIIFGHHCCEEHEHBEDFEFEFD>B#####
@SRR957824.2 2/1
GATCGGAAGAGCACACGTCTGAACTCCAGTCACAGCGCTATCTCGTATGCCGTCTTCTGCTTGA#####CATACTAATATGAGTTCGATATATAAACTGAAAGATCAACATGAGTTCATCAGTGACTATAT
+
???????DDDDDDGGGGGIIIIIIHHIIHHHHHHHHIHHIIHHHHHHIHHGHIIIFGHDHHEHHE#####
@SRR957824.3 3/1
TGGCTAACCCGGATACCACTTCCGGCGAAATGTGCGTATTATCCACAGATTCATCGTTTAACACGAATTTTCAAACGGAACAGCTTATGAGTGCAATCGCGCTGGAATGATCCTCATCGGTACCTCTGCGGCTCCATTTCCAGTGC
+
???A?BBBDDDDDDGGGGGIHHHHHHIIFHHHHHHHEHHIHHHHIIIIIIHHIIIIHHHFHIIHHFHHHE?DHHHHHHHHHFHFGGGGGF=BGGDECCEGGGGGGGGGGC6<>DEEGG?E>D<>EEEEGGEGCGC:(
@SRR957824.4 4/1
AACGCCAGCATTGAGGGACCGCAACCTGGTGAGAACGCTCGGACATCTTAATGGCGACGGCGCTTACTGGAGCGCTTCGCTGGAAGAATCTCAATGTACCGCAGAAAAATCGCGCTTGACGCTGTGGAGCGGGTTATCTGCTTGC
+
,????@?DDDDDDDBEFEFHHHHHHHHHAFGHFFHHHEHACECHHHHHGHHHHHHDHEEAHEHHHEFFFFFEEEEFFEEFFFA;==EEFBEEEBEEEEDD;>?AEC?EFDDEFFFFEECCE?EFCA;;D;AAEA?:A?:?:?
@SRR957824.5 5/1
TGGCGATACGACGGAAGTGATAAGTACCCGTGGTCTCAGCCTGTGCTGTATCACATTTGCCTCGAGATGCGTTCAAGGGGATTGAGCGCCAGATGACCGAAGGGGAATTAACGCGCTTGACGAAACGGCAACTGACGAAATGGGCCAA
+
???AAAAADDEEDDDGGGGGIHHIIIIHHHFGHHIIIIIIHHIIIIIIHHHEHHIIHHHHHHHHHHHHHHHGHDFGG:@EGEGGGGGGGGGEGGGGB'?EGC**8?4'48*?EGGE8C8:CE?*8?
@SRR957824.6 6/1
CTGGCGATTGCTGCTTGGGCGCGGATTACTATCATCATCCGATTGGCGATGTGCTGTTTCATGCCTTGCCGATTTTACTGCGCCAGGGCGGCTGCGGCGAACGCGCCGATGTGGTACTGGTTTGCCACTGAACAAGGCCATGCGC
+
```

```
zcat SRR957824_2.fastq.gz | less
```

**Question 1: Why are there two files for the single sample dataset?**

The next step is to access the quality of our raw reads. We will do this by running FASTQC on our dataset.

```
module load fastqc/0.11.7
```

```
fastqc SRR957824 1.fastq.gz
```

```
fastqc SRR957824 1.fastq.gz SRR957824 2.fastq.gz
```

```
fastqc *.fastq.gz
```

Tip: A tool called MultiQC (<https://multiqc.info/>) can be used to combine the results from several QC outputs into a single report. This may be useful when working with a large number of samples. FASTQC can also be run via a graphical user interface, via the terminal or in R Studio (fastqcr package).

Once the analysis is complete, you will notice that new files with a .html extension will be generated. These files contain the information report and images produced by the FASTQC program. You will need to copy the FASTQC html files back to your local machine in order to view the reports in a web browser.

Open up a new terminal on your local computer (not on the HPC).

Use the following command to copy the html FASTQC files to your local computer. Replace **X** with your username and **/path/to/your/html\_files/** with the actual path to your FASTQC html files on the HPC. You can get the path to your FASTQC html files by typing `pwd` on the terminal while in the directory with the FASTQC html files.

```
scp userX@hpc.ilri.cgiar.org:/path/to/your/html_files/*.html .
```

Once the files have been copied to your local computer, you can open them using the following approaches based on your operating system:

#### MacOS

```
open SRR957824_1_fastqc.html
```

#### Linux

```
firefox SRR957824_1_fastqc.html
```

#### Windows

Navigate to the location where you copied the file to and double-click on the file to open it. If you are using a terminal emulator it should be in the installation directory for the emulator.

This command should open the html report for read1 in your browser. You can use the same approach to open all the FASTQC reports.

# Summary

- ✓ [Basic Statistics](#)
- ✓ [Per base sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ✗ [Per base sequence content](#)
- ! [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ✗ [Overrepresented sequences](#)
- ✓ [Adapter Content](#)


## Summary

- ✓ [Basic Statistics](#)
- ✓ [Per base sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ✗ [Per base sequence content](#)
- ! [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ✗ [Overrepresented sequences](#)
- ✓ [Adapter Content](#)

### Basic Statistics

Measure	Value
Filename	EBI_SRA_SRR957824_File_ftp__ftp_sra_ebi_ac_uk_vol1_fastq_SRR957_SRR957824_SRR957824_1_fastq_gz.gz
File type	Conventional base calls
Encoding	Sanger / Illumina 1.9
Total Sequences	1792335
Sequences flagged as poor quality	0
Sequence length	150
%GC	49

### Per base sequence quality



In a similar way open the quality control report for Read2.

firefox SRR957824\_2\_fastqc.html



# Summary

- ✓ [Basic Statistics](#)
- ✗ [Per base sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ✗ [Per base sequence content](#)
- ✓ [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ! [Overrepresented sequences](#)
- ✓ [Adapter Content](#)


## Basic Statistics

Measure	Value
Filename	EBI SRA_ SRR957824 File_ftp__ftp_sra_ebi_ac_uk_voll_fastq_SRR957_SRR957824_SRR957824_2_fastq_gz.gz
File type	Conventional base calls
Encoding	Sanger / Illumina 1.9
Total Sequences	1792335
Sequences flagged as poor quality	0
Sequence length	150
%GC	50

# Summary

- ✓ [Basic Statistics](#)
- ✗ [Per base sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ✗ [Per base sequence content](#)
- ✓ [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✓ [Sequence Duplication Levels](#)
- ! [Overrepresented sequences](#)
- ✓ [Adapter Content](#)

## ✗ Per base sequence quality



### Task 3: Interpretation of QC report

Compare the quality control reports for the two FASTQ files.

#### Question 2: Comment on the quality of the sequenced reads in file 1 and file 2.

Read through the other reports in the FASTQC output html file and ensure that you understand the information in each report.

#### Question 3: What can be done to improve the quality of the reads?

#### Trimming

Trimming can be done using a tool called Trimmomatic (<https://github.com/usadellab/Trimmomatic>)

The current trimming steps are:

- ILLUMINACLIP: Cut adapter and other illumina-specific sequences from the read.
- SLIDINGWINDOW: Perform a sliding window trimming, cutting once the average quality within the window falls below a threshold.
- LEADING: Cut bases off the start of a read, if below a threshold quality
- TRAILING: Cut bases off the end of a read, if below a threshold quality
- CROP: Cut the read to a specified length
- HEADCROP: Cut the specified number of bases from the start of the read
- MINLEN: Drop the read if it is below a specified length
- TOPHRED33: Convert quality scores to Phred-33
- TOPHRED64: Convert quality scores to Phred-64

It works with FASTQ (using phred + 33 or phred + 64 quality scores, depending on the Illumina pipeline used), either uncompressed or gzip'ed FASTQ. Use of gzip format is determined based on the .gz extension.

Below is an example for using Trimmomatic. Trimmomatic will produce both paired and unpaired datasets from a paired-end set of FASTQ files. The file ILLUMINACLIP:TrueSeq3-PE.fa contains the adapter sequences associated with that sequencing platform.

```
java -jar trimmomatic-0.35.jar PE -phred33 input_forward.fq.gz input_reverse.fq.gz output_forward_paired.fq.gz
output_forward_unpaired.fq.gz output_reverse_paired.fq.gz output_reverse_unpaired.fq.gz ILLUMINACLIP:TruSeq3-
PE.fa:2:30:10 LEADING:3 TRAILING:3 SLIDINGWINDOW:4:15 MINLEN:36
```

This will perform the following:

- Remove adapters (ILLUMINACLIP:TruSeq3-PE.fa:2:30:10)
- Remove leading low quality or N bases (below quality 3) (LEADING:3)
- Remove trailing low quality or N bases (below quality 3) (TRAILING:3)
- Scan the read with a 4-base wide sliding window, cutting when the average quality per base drops below 15 (SLIDINGWINDOW:4:15)
- Drop reads below the 36 bases long (MINLEN:36)

Run the command below in order to trim your fastq files. Ensure that you understand what each of the options in the command means. You may need to load the Trimmomatic module first.

```
trimmomatic PE SRR957824_1.fastq SRR957824_2.fastq  
SRR957824_1_paired.fq.gz SRR957824_1_unpaired.fq.gz  
SRR957824_2_paired.fq.gz SRR957824_2_unpaired.fq.gz LEADING:15  
TRAILING:15 SLIDINGWINDOW:4:20 MINLEN:100
```

**Question 4: Explain what each of the options in the trimmomatic command above specifying. i.e. LEADING, TRAILING, SLIDINGWINDOW, MINLEN.**

If you list your files again you will notice that you now have some new unpaired and paired .fq.gz files. These are your qc'd files that have been trimmed using Trimmomatic in the last step. To see the effect of your trimming on your new files, run FASTQC on the paired output files as you have done previously. Compare the output of the first FASTQC run and the second FASTQC run.

**Question 5: Comment on the quality of your two trimmed read files. Have you improved the quality of your data? How many reads have you lost for each file by trimming the data? Please note any other improvements in the dataset.**